

## Vehicle harness generation 2

- Designed for **children and adults with special needs or challenging behaviour** during transportation.
- Provides **postural support** for people with **physical disabilities** during transportation.

This car harness provides a real feeling of **security** and **comfort** for the users, driver, as well as the other passengers.



### Description:

- **New ergonomic design** with adjustable backrest
- **New emergency evacuation system**, reachable behind the seat or by the trunk
- **Made in France** in neoprene material
- **Adaptable & long lasting** : Can be used with booster seat and has adjustable belts to fit the growing child
- **Universal**: adaptable to all type of seats
- **Safe** : provided with seven lockable buckles and a set of keys
- **Easy to transport** in car, bus, trip, rental car



### Options:

- **Headrest** with **adjustable height** and **front support**
- Set of **2 additional belts** for individual seat installation
- **New emergency evacuation system ESCAPE BOX**
- **Modifications** : **possibility to tailor-made the harness**. Additional cost may apply



**Important information:** The **seat belt remains compulsory** while using the vehicle harness. We recommend you to check the country's regulations with regards to safety car belts prior to using the vehicle harness. The **utilisation** of the vehicle harness remains at the driver's responsibility.

# Installation notice

Support base (A)



Adjustable backrest (B)




Frontal harness (C)



## Installation on standard seat: there is a space between the backrest and the seat

1. Unclip the back support (A+B) from frontal harness (C).
2. Place the back support (A+B) over the backrest of the car seat . If necessary, fold forwards the seat to pass the belts that go around the back of the seat.
3. Clip and tighten the locking belts to the maximum to ensure the harness doesn't float.
4. Once the person is seated comfortably, take the frontal harness part (C) and clip it against the back support (A+ B). Adjust the 7 belts and lock, if necessary.
5. Please return the belt bottoms as per picture



## Installation on individual seat: there is no space between the backrest and the seat

The installation steps are the same as for the standard seat. 2 additional belts are necessary and available as an option for the harness.

1. One additional belt is to fix behind the seat and help the harness back support (A+ B) go around the whole seat, not only the backrest.
2. The second one is to fix between the seat and the backrest. Tighten the belt to the maximum to ensure the harness doesn't float.



**Important information:** The standard emergency evacuation system is reachable behind the seat or by the trunk. You only need to click on the 2 red buttons at the same time and the person can leave the car together with the harness.



We wish you a safe trip!