

Headrest

The headrest provides head support for wheelchair users. When used during transport, the headrest can reduce the risk of whiplash in case of accident.

- **Easily installed** on the wheelchair handles (max. diameter of 25mm).
- **Fully adjustable** for a better comfort (height, width, inclination).
- **Easy care** : covered with PVC (washable)
- Headrest **available in 3 different dimensions**: “Small” (S), “Medium” (M) and “Large” model (L).

(S)



(M)



(L)



Wheelchair straps D rings (DR120)

The “D Ring” is a wheelchair accessory to use in complement to the restraint systems. **It allows the restraint to be secured into place.** The “D Ring” accessory is tested to **ISO10542 requirements.**



Secure straps for transport of an empty folded wheelchairs (RS)

Easy to fasten and install, the wheelchair strap enables to secure **empty folded wheelchairs** and other products **during transport.**

